

PRODUCT BULLETIN #3231

To: Our Valued Customers

From: Michael Verguldi, Sr. Product Manager

Date: October 16, 2020

Subject: **PRS-12 Discontinued**



INTRODUCTION

The PRS-12 has been discontinued effective immediately and is no longer available. It will be replaced with the PRS-13 System.

BACKGROUND

The PRS-12 System, part number 3028096, which was made available to the market in August of 2012, has been discontinued due to the manufacturer obsoleting some of the key components. As a result, we have released a new Particle Reduction System, the PRS-13.

REPLACEMENT SYSTEM

The PRS-13 operates the same as the PRS-12 and has many of the same features and benefits including performance, design, operation, and overall size. The key difference is the PRS-13 has a different ultra-filter (UF) housing with end caps. A new brochure will be available for the PRS-13 that details product description, specifications and technical data.

TRANSITION AND SUPPORT

The PRS-12 System will continue to operate as intended and Mar Cor will continue to support this system, but there are some exceptions. Many of the PRS-12 spare parts are interchangeable with the new PRS-13 system, but there are a few that are no longer available. If the housing or end caps become damaged or require replacement, the PRS-13 system will need to be purchased instead since the housing and/or associated housing fittings are no longer available. The ultra-

filter cartridge itself is still available and is designed to work in either system. See the chart below for more information:

Part Number	Description
Obsolete Parts	
3028065	HOUSING, ULTRA-FILTER CARTRIDGE, PRS-12
3028100	END CAP KIT, PRS-12
3028101	LOCKING RING, END CAP, PRS-12
Interchangeable Parts	
3028093	ULTRA-FILTER CARTRIDGE, REPLACEMENT KIT, PRS-12/13
3028097	SOLENOID VALVE, SYSTEM FLUSH, PRS-12/13
3028098	TRANSFORMER, 24V, PLUG-IN ADAPTER, PRS-12/13
3028099	PULSE CONTROLLER, PRS-12/13
3028064	PRESSURE TANK, PRS-12/13

We thank you for your understanding and value your business. If you have any questions, please contact your local Sales Representative or Product Management at 800-633-3080 or visit our website at www.mcpur.com.