

ELGA® LabWater Products

Mar Cor Purification is proud to be a distributor of ELGA LabWater products. ELGA LabWater, a trusted name in water purification systems, is now teamed up with Mar Cor Purification's extensive sales and service network in North America.

ELGA LabWater systems meet the needs of research, science, clinical and healthcare environments with pure and ultra pure water purification systems from just a few liters up to several thousands of liters of water per day. Described below is a brief overview of the portfolio of products available.



PURELAB®

Purified Water for Laboratory Research and Testing Applications

From simple routine washing and rinsing through to ultra pure water being used for the most critical applications such as ICPMS and ultra trace analysis, the PURELAB range is designed to offer water purity to match your specific laboratory application. The PURELAB product line offers a range of products from Type I to Type III and can produce flow rates of 1 to 1000 liters per day.

Popular PURELAB Products

PURELAB Classic

- The best choice for budget minded laboratories needing top quality water
- Consistent supply of 18.2 MΩ-cm water, at flow rates of up to 2 liters per minute

PURELAB 3000

- High volume delivery of pure water for general laboratory use
- Flexible configuration provides flow rates from 60 to 120 liter per hour of Type II Water

PURELAB Ultra

- Choice of three ultra pure water systems for chemical analysis applications
- Purity beyond 18.2 MΩ-cm with unique PureSure system

PURELAB flex

- Uniquely designed dispensing system delivers high quality water at up to 2 liter per minute
- Point of use display indicates water purity as you dispense





MEDICA®

Purified Water for Clinical Diagnostics

MEDICA products are designed to provide water to clinical analyzers specified to Clinical Laboratory Reagent Water (CLRW) standard (formally Clinical Laboratory Standards Institute Type 1). Whether you have access to tap or partially purified water, or need to feed analyzers on a small or large scale, ELGA has the system for you.

MEDICA-R 7/15

- Delivers up to 1.8 liters per minute of CLRW standard water (formally CLSI Type 1 water)
- Provides compliant purified water suitable for a dissolved oxygen specification

MEDICA Pro

- Delivers up to 4 liters per minute of CLRW standard water
- Compatible with single and multiple clinical diagnostic analyzers which require up to 4 liters per minute of CLRW standard water

MEDICA-R 200

- Delivers up to 21 liters per minute of CLRW standard water
- Can feed multiple analyzers on a distribution loop and across several laboratories
- A bypass loop for continual supply of compliant water in emergency situations
- Compact footprint due to its integrated wrap around 350 liter reservoir
- Compatible with multiple and large clinical diagnostic analyzers which require up to 21 liters per minute of CLRW standard water

CENTRA®

Pure Water Systems for Centralized Distribution of Laboratory Water

CENTRA products are single compact units that purify, store and distribute water from a single central unit. These pre-designed standard products are ideal for projects requiring a higher quantity of water packaged in a cost effective design.



CENTRA R 60/120

- Complete water purification, storage, control and distribution systems
- Delivers up to 60-120 liters per hour (0.25-0.5 gpm) of makeup water
- Delivers up to 2.5 gpm of distributed water
- Compact design means that the system can be placed close to work areas

CENTRA S 200, R 200

- Centralized purification and distribution systems capable of supplying pure water to one or more laboratories on a loop system
- Capable of delivering 200 liters of pure water per hour (0.9 gpm) of makeup water and 4.7 gpm of distributed water from integral tank

ELGA is the global laboratory water brand name of Veolia Water Solutions & Technologies. ELGA®, PURELAB®, MEDICA®, CENTRA®, ADEPT® Labpure® and PureSure® are registered trademarks of VWS (UK) Ltd.



Mar Cor Purification
4450 Township Line Road
Skipack, PA 19474-1429
Tel: (484) 991-0220
Toll Free: (800) 633-3080
Fax: (484) 991-0230

Mar Cor Purification
14550 28th Avenue North
Plymouth, MN 55447
Tel: (484) 991-0220
Toll Free: (800) 633-3080
Fax: (763) 210-3868

Mar Cor Purification
1119 Paulsun Street
San Antonio, TX 78219
Tel: (210) 227-3601
Toll Free: (800) 268-5035
Fax: (210) 227-0735

Mar Cor Purification
160 Stedman Street
Lowell, MA 01851
Tel: (978) 453-9600
Toll Free: (800) 633-3080
Fax: (978) 453-1223

Mar Cor Purification
6351 Orangethorpe Ave.
Buena Park, CA 90620
Tel: (714) 736-9990
Toll Free: (800) 633-3080
Fax: (714) 736-9402