

Mar Cor Purification Life Science Products



Mid-Atlantic Service Centers



Mar Cor Purification Life Science Products



Corporate Profile

Mar Cor Purification (MCP) is owned by Cantel Medical Corporation, a New York Stock Exchange listed company (CMN). MCP is a leading provider of complete water treatment systems and services for the Life Science market. For over 40 years and over 1000 installations, we have been providing high purity water solutions to these industry leading sectors. The Philadelphia Service Center is one of 27 MCP offices and one of 6 regeneration plants across the US and Canada. The Philadelphia Service Center contains a SDI regeneration plant that is registered with the FDA as a medical device manufacturer, plus features a qualified team of Mar Cor Purification Field Service Technicians.

Mar Cor Purification Water Treatment Capabilities

- Reverse Osmosis (RO)
- Water Softeners
- Carbon Filters
- Media Filtration
- Electrodeionization (EDI)
- Service Deionization (SDI)
- Ultrafiltration
- Ultraviolet Purifiers
- Depth/Pleated Filtration
- Ozone Disinfection Systems
- Cleaners & Sterilants
- PLC/HMI Control Systems



Pretreatment Systems

Pretreatment for Lab Research and Medical Applications



BioPure LSX

High Efficiency USP High Purity Water System

Products & Services

High Purity Water Systems

Mar Cor Purification (MCP) is an industry leader in high purity water systems for Medical, Pharmaceutical, Semiconductor, Commercial & Industrial, and Lab Water applications. MCP can provide a quality water purification solution for any situation or application.

- Single-source provider of installed water purification systems
- Design & Engineering Services
- Standard and custom modular systems, including skid mounted
- Integrated PLC based control systems
- Superior designs proven to minimize bacteria growth
- "Green" options for water recovery and energy efficiency



Typical Applications

- Pharmaceutical Preparations
- Biological Engineering
- Central Laboratory Water
- Power Industry
- Humidification
- Rinsing of Electronic Components
- Glassware Rinsing
- Preparation of Process Chemicals

Standard Features

- Skid Mounted
- Integral Power Supply
- Operational Instruments
- Epoxy Coated Steel Frame
- Pressure Gauges
- Flow Meters
- Prewired/Preplumbed
- Factory Tested

Optional Items

- Hot Water Sanitizable
- Various Piping Materials
- Stainless Steel Frame
- Brine Feeder
- Ultraviolet Sterilizer
- NIST Traceable Instrumentation
- High Water Recovery Design
- System Validation



Service Deionization (SDI)

Mar Cor Purification offers portable exchange activated carbon and deionization services. All portable exchange medias are processed in our FDA registered and compliant facilities.

- No capital investment required
- Monthly, Quarterly or Semi-Annual Services available
- Premium grade carbon and ion exchange resin medias (cation, anion or mixed bed)
- Flow rates from 0.25-100 gpm. Additional flow rates available with parallel units
- FDA Registered, Compliant and Consistent



Skid Mounted Purification Equipment

- Pre-Engineered Models
- Complete Skid Mounting Options
- Compact Footprint
- Pre-Piped and Wired Controllers
- Reduce On Site Installation Times
- Performance Testing Options Available

Mar Cor Purification Support Services

Mar Cor Purification (MCP) is one of the largest and most experienced water purification companies in North America. We have been servicing water treatment equipment for over 40 years and have an experienced team of technicians qualified to work on your water systems based out of our Skippack, PA office. Key services offered include:

Engineering and Design Services

Working with our clients, in a collaborative engineering and design effort, ensures the delivery of equipment properly designed to meet your needs. MCP provides fully integrated services so that the system you purchased meets your requirements.

After Care Services

Once installed, MCP preventative maintenance contracts and services safeguard the continued production of high purity water. MCP is available 24/7 to provide you emergency and scheduled support services to keep your system operating at peak efficiency including:

- Periodic service by qualified technicians to ensure consistent performance and protection of your capital investment
- Analytical testing by independent laboratories to stay compliant
- Periodic system cleaning and sanitization to maintain high quality standards
- Regular equipment calibrations to NIST standards

Equipment Installation

- System Start-Up
- In-Depth Operator Training
- Storage & Distribution Piping Systems

Maintenance & Repair

- Diagnostics and Repair Services
- System Performance Audits
- Preventive Maintenance

Engineering Services

- Validation Documentation
- Engineering Documentation
- Project Management



Local Contact Information for Philadelphia Service Center:

Mar Cor Purification
4450 Township Line Rd
Skippack, PA. 19474
Phone: 484-991-0220



Mar Cor Purification 4450 Township Line Road Skippack, PA 19474-1429 Tel: (484) 991-0220 Toll Free: (800) 346-0365 Fax: (484) 991-0230	Mar Cor Purification 14550 28th Avenue North Plymouth, MN 55447 Tel: (484) 991-0220 Toll Free: (800) 633-3080 Fax: (763) 210-3868	Mar Cor Purification 3250 Harvester Road - Unit 6 Burlington, ON L7N 3W9 Tel: (905) 639-7025 Toll Free: (800) 268-5035 Fax: (905) 639-0425	Mar Cor Purification Sourethweg 11 6422 PC Heerlen The Netherlands Tel: (+31) 45 5471 471 Fax: (+31) 45 5429 695	Mar Cor Purification 1A International Business Park, #05-01 Singapore 609933 Tel: (+65) 6227 9698 Fax: (+65) 6225 6848
---	--	---	---	---