

MAR COR PURIFICATION LAUNCHES NEW USP PURIFIED WATER SYSTEM- BIOPURE LSX

... New Design Features Energy Conservation and Fully Automated Hot Water Sanitization

Plymouth, MN - January 22, 2013 - Mar Cor Purification has launched a new USP Standard Water System design, the BioPure LSX. The BioPure LSX system is a single skid, single pass/RO/EDI USP Purified water system, used in pharmaceutical, cosmetics, medical device and other high purity/bacteria controlled applications.

Ben Rocznik, VP Sales, commented, "The BioPure LSX is the latest water system designed to address the needs of the Biotech and Pharma water markets. The previously launched US Pure and Vpure designs were each designed to meet the needs of a specific market segment. The BioPure LSX is designed to install quickly, incorporate energy saving technology, automated hot water sanitization of the system, and a high level of automation to adjust operating parameters to optimize the operation. All of these features are offered at a very competitive price point. The standard systems are available in flowrates up to 20 gpm. Mar Cor's leadership position in high purity water enables us to introduce innovative designs that meet the current needs of the marketplace."

About Mar Cor Purification

A subsidiary of Cantel Medical Corp., Mar Cor Purification is dedicated to providing innovative solutions through filtration, water, and disinfection technologies to the medical, life science and industrial marketplaces. Mar Cor Purification has service offices in 19 cities in the U.S. and Canada, with 7 resin regeneration plants strategically located in Atlanta, Boston, Chicago, Philadelphia, San Antonio, Montreal and Toronto. For more information please visit www.mcpur.com.

About Cantel Medical Corp.

Cantel Medical Corp. (NYSE:CMN) is a leading provider of infection prevention and control products in the healthcare market. Our products include water purification equipment, sterilants, disinfectants and cleaners, specialized medical device reprocessing systems for endoscopy and renal dialysis, disposable infection control products primarily for dental and GI endoscopy markets, dialysate concentrates and other dialysis supplies, hollow fiber membrane filtration and separation products for medical and non-medical applications, and specialty packaging for infectious and biological specimens. We also provide technical maintenance for our products and offer compliance training services for the transport of infectious and biological specimens.

This press release contains forward-looking statements within the meaning of the Private Securities Litigation Reform Act of 1995. These statements involve a number of risks and uncertainties, including, without limitation, the risks detailed in Cantel's filings and reports with the Securities and Exchange Commission. Such forward-looking statements are only predictions, and actual events or results may differ materially from those projected or anticipated.

SOURCE: Mar Cor Purification

Contact information:

Christopher Fournier

cfournier@mcpur.com

Vice President Marketing

Mar Cor Purification

Phone: 978-453-9600