

Millenium HX™ Water Purification System Check Valve Preventative Maintenance

To reduce machine downtime, an annual preventative maintenance (PM) task has been added to the Millenium HX RO machine as of December of 2019.

The PM kit (P/N 3035238) is specifically designed to replace three check valves on the Millenium HX: Check Valve Five (CK5), Check Valve Seven (CK7), and Feed Check Valve (FCV).

The Check Valve PM kit is intended to simplify this recommended annual work as it includes all the necessary check valves and written instructions on proper replacement in one easy-to-order kit.

Please contact Technical Support at 800-633-3080 if you have any questions.

Millenium HX™ is a trademark of Mar Cor Purification, Inc.