

Millenium RO Water and Power Consumption

This table represents the typical operational flow and power consumption of the Millenium RO at different operation conditions. These values assume that the machine is operating within the parameters of the O&M manual.

Water and Power Characteristics of the Millenium RO	
Operational flow @ 800 mL/min (13°C)	1.0 - 1.1 GPM 3.8 - 4.2 LPM
Operational flow @ 500 mL/min (13°C)	0.80 - 1.0 GPM 3.0 - 3.8 LPM
Operational flow without patient draw (13°C)	0.53 – 0.75 GPM 2.0 – 2.8 LPM
Operation flow without dialysis machine (13°C)	1.3 - 1.5 GPM 4.9 - 5.7 LPM
Typical Electrical Consumption	560 - 630 Watt