

MILLENIUM REPAIR/CHECK TEST DATA

Model Number: _____ **Service Date:** _____
Part Number: _____ **Technician:** _____
Serial Number: _____

Description of Initial Problem:

Description of Service Performed:

Was Problem Corrected ? Yes No

Millenium Check Test Data	Actual	Specification
Feed Water Temperature	_____	77°F (25°C)
Pre-Filter Pressure (psi)	_____	≥30 PSI
Post Filter Pressure (psi)	_____	≤ 100 PSI
Pump Pressure (psi)	_____	190 PSI (±10 PSI)
Product Pressure (psi) – product water hose running to drain	_____	1 PSI
Prod Pressure (psi) – product water hose disconnected	_____	30 PSI
% Rejection	_____	
Product Water TDS (machine)	_____	
Product Water TDS (manual)	_____	
Product Water Flow *	_____	0.52 gpm (1.97 lpm)
Waste Water Flow	_____	0.84 gpm (3.2 lpm) +/-10%

* Prod water flow rate specific to machine PN 1232429. Product water flow rates for the earlier machines (PN 14714) 0.34 gpm (1.31 lpm).