
Minncare[®] Cold Sterilant and Minncare[®] Test Strip Expiration Dates

Customers have asked about the expiration date once Minncare[®] Cold Sterilant and Minncare[®] Test Strip containers have been opened.

Minncare[®] Cold Sterilant

The expiration date of Minncare Cold Sterilant remains unchanged from the expiration date printed on the label, as long as the Minncare Cold Sterilant is stored in the original container, the container is closed using the original vented container cap, the container is stored in a cool, dry place (below 75°F), and is not exposed to direct sunlight.

Product that has been dispensed from the original container and is not used should not be returned to the container, but should be discarded in accordance with all federal, state and local regulations.

Minncare[®] 1% and Residual Test Strips

The expiration date listed on the Minncare 1% and Minncare Residual Test Strips list the month and year of expiration, this expiration date is an “inclusive” date.

As long as the test strips are properly stored and handled, they are valid and can be used through the last day of the month listed on the container, even after opening.