

PRODUCT BULLETIN #3230

To: Our Valued Customers

From: Michael Verguldi, Sr. Product Manager - Medical

Date: August 24, 2020

Subject: **CWP 100, WRO H/S Operator's Panel**



INTRODUCTION

We qualified a new operator's panel for the water purification unit, CWP 100, model WRO H/S.

BACKGROUND

The operator's panel, also referred to as a Human Machine Interface (HMI) controller, used in the CWP 100 models WRO H/S has been discontinued due to the manufacturer's decision to no longer offer the E1032 line of controllers. As a result, we have released a new operator's panel for all CWP 100 units.

REPLACEMENT OPERATOR'S PANEL

The replacement operator's panel is a full-color touch screen and is designed to have a similar look and feel as to the existing E1032 controller. Additional information is found in Tech Note 312.



Figure 1: X2 Base 7 (new)



Figure 2: E1032 (existing)

TRANSITION

The current CWP 100 E1032 operator's panel will be discontinued as of October 5, 2020 and will be replaced with the X2 base 7 HMI. All CWP 100 WRO H/S models built on or after that date will have the new X2 base 7 HMI. For installed CWP 100 units, Mar Cor will maintain an inventory of spare E1032 controllers until supplies last. Once those controllers are no longer available, designated replacement kits have been created to either change a E1032 that requires service or to simply upgrade to the new X2 base 7 operator's panel. The chart below is designed to assist with obtaining the proper kit part number for either objective:

Models H/S	Part Number	Description
101-102 H	3035783	KIT, HMI, UPGRADE, X2 BASE, CWP, 101-102H
104-106 H	3035785	KIT, HMI, UPGRADE, X2 BASE, CWP, 104-106H
101-102 S	3035782	KIT, HMI, UPGRADE, X2 BASE, CWP, 101-102S
104-106 S	3035784	KIT, HMI, UPGRADE, X2 BASE, CWP, 104-106S

We thank you for your understanding and value your business. If you have any questions, please contact your local Sales Representative or Product Management at 800-633-3080 or visit our website at www.mcpur.com.